

# WASTE STATION

## Stand-alone bedding disposal station



### Introduction

The requirement for a bedding disposal station in your facility is clear and simple. It is a neat way of ensuring a clean working environment. Making life easier for animal technicians to dispose of waste, quickly and efficiently.

### Active Airflow Technology

The system draws air from the room into the work area. Just above the work surface is where the air has a maximum flow, with a further option to extract air from the bin store below. Two levels of pre-filtration come as standard, situated behind a protective baffle. A cleanable fine mesh prevents ingestion of larger particles; then washable prefilters contain smaller dust particles. The final stage is an F9 Microfilter, limiting your exposure to allergens, odours and waste while emptying dirty cages. Further to this, the system continually monitors airflow to alert the user if the filters are getting blocked, or whether a filter is damaged.

Dimensions			
Dimension Description	Inches	mm	
Width	45.43	1154	
Height	77.17	1960	
Depth	29.33	745	
Largest Cage (Height)	13.78	350	
Bin Capacity	Cubic Inch	CC	Litres
Full lined bin	7323	120000	120

Weight		
Unpacked Weight	lbs	kgs
	419	190

Noise	
Noise	dB
Fans only	57



ZOONLAB GmbH  
Hermannstraße 6  
44579 Castrop-Rauxel  
Deutschland



T: +49 (0) 23 05 - 9 73 04 - 0  
F: +49 (0) 23 05 - 9 73 04 - 44  
E: info@zoonlab.de  
www.zoonlab.de

## Improved Efficiency

The station's ergonomic design lends itself to being extremely user-friendly and easy to clean. All components are detachable and can be placed in a cage washer- saving you time and money on having to manually deep clean soiled components. For everyday handling, each unit comes as standard with a magnetic scraper to clean down the unit after each use.

## Dual Filtration

Sychems' stand alone bedding disposal station enables a cleaner working environment for all. Utilising Active Airflow Technology, a filtered laminar airflow creates a downflow of H14 HEPA

## Unrivalled Cleanability

Simply, the entire work platform can be lifted out for manual cleaning. The clear plastic sides are easy to remove, with both the chute and worksurface able to lift out. Making it easier than ever for you to clean after use.

## Easy On-Site Set Up

Powered via a standard mains socket, the standalone waste disposal unit comes standard with a single chute and tabletop. Providing you with a spacious working area, compatible for working with all cages both small and large.

The station is lightweight and can be positioned into the desired operation seamlessly. The unit will fit through any standard door and can be secured safely with lockable front wheels.

A neat and tidy display screen is positioned in the right-hand corner of the unit, displaying an information status screen. Visual identification through a colour-coded emoji-styled set of indicators. Green/Amber/Red symbols highlight whether the system's filtration system is operating efficiently.

## Cage Compatibility

The unit's spacious surface area provides you ample room to work with cages of all shapes, sizes and weights. You can rest assured knowing that SychemWASTE will provide your facility with the ideal solution for maintaining a clean and tidy working environment.

## Standard Features

- One tabletop, chute, and bin
- Removeable knockout bar
- Three stage filtration
- Two-year warranty
- Magnetic scraper

## Accessories

- Second chute
- Sychem bin(s)
- Side handles
- Exhaust thimble
- Magnetic bumpers

