

CHANGESTATION

Dual-sided animal transfer station

NEW



With quiet operation and full night mode control, SychemCHANGE minimises the impact on the animal's sleep cycle as well as keeping stress levels at a minimum. Height adjusting legs allow the SychemCHANGE unit to be moved to a comfortable operating height for all users.

Night mode is a user-set feature that automatically sets the LED lighting to red and changes the touchscreen interface to a red graphical overlay at a specified time. This allows the machine to be operated during the animal's sleep cycle, without disrupting the animals. SychemCHANGE has been designed with an intuitive user interface, allowing for easy configuration and control over all parameters, while keeping relevant information a glance away.


What is an animal transfer station

Keeping animals happy, healthy and stress-free means moving them to a fresh, clean cage regularly. In order to keep both the animals and the humans safe, the animals are moved to the clean cage under a filtered laminar airflow, keeping all animal dander and allergens away from the technician and keeping any foreign bacteria and dander away from the animals. SychemChange does this by creating a downflow of H14 HEPA filtered air down on the work surface, this air is then sucked down through the work surface, filtered again and then exhausted out the bottom. This method of using dual filters reduces the build-up of heat and odours, as well as elongating the filters life.

Dimensions		
Dimension Description	Inches	mm
Width	63.46	1612
Height	64.52	1639
Raised Height	82	2089
Depth	30.39	772
Largest Cage (Height)	15.74	400

Weight		
Weight	lbs	kgs
Unpacked Weight	580	263

Noise	
Noise	db
Fans only	55



ZOONLAB GmbH AHE
Hermannstraße 6
44579 Castrop-Rauxel
Deutschland



T: +49 (0) 23 05 - 9 73 04 - 0
F: +49 (0) 23 05 - 9 73 04 - 44
E: info@zoonlab.de
www.zoonlab.de

Height Adjustment

Height adjusting legs allow the SychemCHANGE unit to be moved to a comfortable operating height for all users.

Intuitive UI

SychemCHANGE has been designed with an intuitive user interface, allowing for easy configuration and control over all parameters, while keeping relevant information a glance away.

Night Mode

Night mode is a user-set feature that automatically sets the LED lighting to red and changes the touchscreen interface to a red graphical overlay at a specified time. This allows the machine to be operated during the animal's sleep cycle, without disrupting the animals.

Optional Accessories

CHANGEstation has several features as added extras, these include:

- Side-mounted shelves (of various designs)
- Food hopper
- Automatic foaming hand sanitiser holder
- Repositionable locking arm with the following attachments:
 - Magnifying glass
 - Phone holder
 - Clamp
- Rear cover for single-sided operation.

Well-Being

With quiet operation and full night mode control, SychemCHANGE minimises the impact on the animal's sleep cycle as well as keeping stress levels at a minimum.

Easy Access

With its simple design and modular internals, SychemCHANGE can be quickly and efficiently maintained, and should any components fail; these can be quickly replaced, minimising machine downtime.

