

TRANSPORT SYSTEM FOR WATER BOTTLES

TRANSPORT TROLLEY AND BOTTLE WASHING BASKETS

The Z0ONLAB transport system for bottle washing baskets makes all of the work processes associated with drinking bottles easier on a day to day basis.

The transport trolley and the matching washing basket ensure an ergonomic and efficient water bottle management.





TRANSPORT TROLLEYS

FOR BOTTLE WASHING BASKETS, MADE OF STAINLESS STEEL

ZOONLAB transport trolleys for washing baskets make the day-to-day work with water bottles easier. The loaded washing baskets remain in the transport trolley for the entire procedure during the loading and unloading process. Plastic rails (option) simplify the removal of the baskets and prevent bothersome grinding noises. Round tubes and autoclavable easy-running castors ensure excellent transport comfort. The rack made of stainless steel is glass-bead blasted and autoclavable. The front two of the four wheels can be locked so that the user is able to switch between precise mobility and a firm standing. The standard model boasts capacity for 10 bottle washing baskets but transport trolleys for 12, 15 or 18 baskets can also be produced on request.



Type	Description	Dimensions in mm (W x D x H)	Order no.
TW-FSPÜ-PA	holds 10 bottle washing baskets, each 250 or 270 mm high, incl. guide rails with polyamide attachment	610 x 640 x 1,631	5113705-8
TW-FSPÜ	holds 10 bottle washing baskets, each 250 or 270 mm high, excl. guide rails with polyamide attachment	610 x 640 x 1,631	5113705-9
Option			
TW cover hood	Cover hood to protect the transport trolley and its contents, complete with 2 zips, autoclavable	Suitable for the aforementioned transport	3119947

BOTTLE WASHING BASKETS

MADE OF STAINLESS STEEL

ZOONLAB bottle washing baskets are designed for all conventional-sized water bottles. They are made of stainless steel, are electrolytically polished and autoclavable. All of the washing basket versions can be stacked amongst each other. The bottles are fixed into place via a detachable cover which also ensures that they remain secure during the cleaning process.



Type	Description	Dimensions in mm (W x D x H)	Order no.
FSPÜ 25-30 / 250	Washing basket which holds 18 F 25 and/or F 30 type bottles	500 x 265 x 250	5122099
FSPÜ 40-75 / 250	Washing basket which holds 18 F 40, F 50, F 53 and/or F 75 type bottles	500 x 265 x 250	5122100
FSPÜ 25-30 / 270	Washing basket which holds 18 F 25 and/or F 30 type bottles	500 x 265 x 270	5122096
FSPÜ 40-75 / 270	Washing basket which holds 18 F 40, F 50, F 53 and/or F 75 type bottles	500 x 265 x 270	5122097
Option			
FSPÜ cover	Suitable for type FSPÜ 25-30 and FSPÜ 40-75 washing baskets	495 x 240, handle 70 mm high	5122098

Subject to technical modifications.