

## ENRICHMENT

### OCCUPATION FOR LARGER ANIMALS

---

Our polycarbonate enrichment is ideal for scientific animal housing. The material is tough, autoclavable and can withstand the forces of larger animals. Foraging games in particular keep animals busy and give them a sense of achievement.





# FEEDING GAME – ENRICHMENT

## FOR LARGE ANIMALS




### Feed me

- For occupation
- Can be filled with dry food
- Autoclavable up to 134°C
- Non toxic
- Also suitable for large animals such as pigs, sheep or horses, dogs, cats or rabbits
- Dimensions (L x W x H): 270 x 110 x 110 mm
- Material: Polycarbonate
- 1 piece / Packing unit
- Order no. 3084300

### ZOONLAB developed a posture for mounting on enclosure parts

A bracket has been developed to allow the feed roller to be installed in a fixed position on enclosure systems.



### Feed me with bracket

- For occupation
- Can be filled with dry food
- Autoclavable up to 134°C
- Non toxic
- Also suitable for large animals such as pigs and sheep or horses
- Material of the bracket: stainless steel
- 1 piece / Packing unit
- Order no. on request

