


STAINLESS STEEL ENRICHMENT FOR RODENTS

Our stainless steel enrichment is ideally suited for scientific animal housing. This sustainable alternative is readily accepted by the animals. Cleaning, washing and sterilisation processes have almost no impact on product life. A clean, sustainable alternative to our plastic enrichment range.






STAINLESS STEEL – ENRICHMENT FOR RODENTS

MOUSE AND RAT ARCHES

stainless steel AISI 304



Mouse arches reduce aggression between animals and allow them to exercise their need to climb.

The round shape of the arches has been proven in scientific studies and is readily accepted by the animals. Natural behaviour is encouraged.

Type	Dimensions (L x W x H) in mm	Order no.
small	142 x 90 x 60	5083804
medium	160 x 95 x 70	5083801
large	166 x 90 x 94	5083805

MOUSE RAMP



The stainless steel mouse ramp can be inserted into cages as an enrichment. The ramp can ensure accessibility to food and water for animals. Additionally, it serves as a basement enrichment that encourages exploration.

Type	Order no.
For cages Type 2 L	5084130

STAINLESS STEEL – ENRICHMENT

FOR RODENTS

MOUSE-CLIMBER



Mouse Climber is easy to fix on the lid of the cage and creates another climbing possibility for the animals. This makes it easier for the animals to reach otherwise inaccessible places.

Type	Dimensions (L x W x H) in mm	Order no.
Mouse-Climber	142 x 70	5063299