




[COSY MAT]



Originally developed by a working group at Charité, ZOONLAB took up the concept of the injection tunnel and developed a stainless steel version.

Durability and hygiene were key considerations, so all components are washable and autoclavable.

The **injection tunnel, including a cozy mat**, is available under **order number 3131085**.

Additional **replacement mats** can be ordered under **order number 3131088**.

User experience shows that mice like to burrow into the cozy mat, making injections stress-free.

The injection tunnel offers:

- **Protection and a place of retreat for the animals**
- **A cosy zone for injections**
- **High stability**

As always, ZOONLAB focuses on animal welfare.

Dimensions: 140 x 100 x 111 mm (B x L x H) Ø 84,9 inside
made of stainless steel

order number: 3131085

Optional accessories

additional mat

order number: 3131088