


plug-in rat house



ZOONLAB, as a solution finder, has responded to its customers' desire to offer more enrichment products for the humane keeping of rats.

ZOONLAB took this as an opportunity to develop a suitable rat house in a plug-in form. It is a durable alternative to polycarbonate or disposable cardboard products.

Thanks to its simple plug-in connections, it can be assembled quickly and easily and can also be cleaned and autoclaved in the smallest of spaces.

The rat house offers:

- Protection and a place of retreat for the animals
- Two exits for stress-free interaction between the animals
- A cosy area

As is so often the case, sustainability and animal welfare were at the forefront of this development.

Dimensions: 150x150x100mm; Trespa 3mm

Order number: 5083751