

POLYCARBONATE ENRICHMENT FOR RODENTS

Our polycarbonate enrichment is ideally suited for scientific animal husbandry. Inside the polycarbonate enrichment, the animals feel protected and safe whilst also being clearly visible to the users.





POLYCARBONATE ENRICHMENT

FOR RODENTS



Mouse shelter

- Colour: Red transparent
- Autoclavable, non-toxic
- Stackable
- Provides a feeling of safety and protection
- Dimensions: 100x90x55 mm (L x W x H)
- Order no. 3084010



Mouse igloo with running plate

- Mouse-igloo plate: red-transparent
- Running plate colour: Amber
- Autoclavable, non-toxic
- The running plate promotes the instinctive running behaviour and increases activity
- The products can also be purchased individually.
- Mouse-Igloo Order no. 3084011
- Running plate Order no. 3084012



Tunnel

- Material: Polycarbonate - 4 mm
- Autoclavable
- Promotes animal activity, reduces boredom and aggressive behaviour
- Mouse tunnel dimensions: 100 x Ø 50 mm
- Order no. 3084040
- Rat tunnel dimensions: 155 x Ø 75 mm
- Order no. 3084014




Crawl Ball

- Colour: Red transparent
- Autoclavable, non-toxic
- Provides the animal with a good opportunity to hide
- Order no. 3084042



Guinea pig house

- Colour: Red transparent
- Autoclavable, non-toxic
- Option to withdraw
- Dimensions: 205 x 157 x 115 mm (L x W x H)
- Order no. 3084051



Clip on handling tube

- Colour: transparent or red
- easy to fix on the lid
- Autoclavable, non-toxic
- Option to withdraw
- Length 110 mm
- Material: Polycarbonate
- 1 piece per Packing unit
- Order no. 3084193 transparent
- Order no. 3084247 red

POLYCARBONATE ENRICHMENT

FOR RODENTS



Tunnel Clip

- Easy to fix on the lid
- Autoclavable, non-toxic
- Material: Plastic
- 100 pieces / packing unit
- Order no. small 3084157
- Order no. large 3084159



Mouse tent

- Colour: Red - transparent
- Autoclavable, non-toxic
- Provides the animal with a good opportunity to hide
- Dimensions (L x W x H): 130 x 130 x 38 mm
- Material: Polycarbonate
- 1 piece / packing unit
- Order no. 3084086



Nombraero

- Colour: Red
- Easy to fix on the lid
- Autoclavable, non-toxic
- For placing food
- Material: Polycarbonate
- 1 piece / packing unit
- Order no. 3084137



C-Saw

- For Occupation
- Autoclavable, non-toxic
- offers protection and promotes activity
- Dimensions: (L x W x H)
- Material: Polycarbonate
- 1 piece / packing unit
- Order no. 3084195



Clear Handling Tubes

Grabbing mice by the tail should always be avoided. The use of handling tubes makes the transfer of mice stress-free for people and animals.

- Autoclavable, non-toxic
- Reusable
- Mice quickly get used to the tunnels
- Material: Polycarbonate
- 1 piece / packing unit
- Order no. 100 mm length 3084104
- Order no. 130 mm length 3084094




POLYCARBONATE ENRICHMENT

FOR RODENTS




Mouse Mezzanine

- Developed by users
- Autoclavable, non-toxic
- Reusable
- Material: Polycarbonate
- 1 piece / packing unit
- Order no. 3084130




Rodent Retreat

- Autoclavable, non-toxic
- Reusable
- Offers protection and promotes activity
- Colour: red and blue
- Dimensions (L x W x H): 95 x 95 x 152.5 mm
- Material: Polycarbonate
- 1 piece / packing unit
- Order no. red: 3084063
- Order no. blue: 3084219



Fat Rat Hut

- Autoclavable, non-toxic
- Reusable
- offers protection and promotes activity
- Dimensions (L x W x H): 150 x 165 x 85 mm
- Material: Polycarbonate
- 1 piece / packing unit
- Order no. red: 3084059



Mouse Wedge

- Autoclavable, non-toxic
- Reusable
- offers protection and promotes activity
- Dimensions (L x W): 105 x 105 x 150 mm (H) 38 mm
- Material: Polycarbonate
- 1 piece / packing unit
- Order no. 3084087



Mouse/Rat Hang

- Stainless steel
- Reusable
- An alternative for hanging handling tubes on the lid
- Material: stainless steel
- 1 piece / packing unit
- Order no. 4133005 (Rat)
- Dimensions (L x H): 165 x 90 mm
- Order no. 4133006 (Mouse)
- Dimensions (L x H): 140 x 37 mm