

ENRICHMENT

OCCUPATION FOR LARGER ANIMALS

Our polycarbonate enrichment is ideal for scientific animal housing. The material is tough, autoclavable and can withstand the forces of larger animals. Foraging games in particular keep animals busy and give them a sense of achievement.





FEEDING GAME – ENRICHMENT

FOR LARGE ANIMALS





Feed me

- For occupation
- Can be filled with dry food
- Autoclavable up to 134°C
- Non toxic
- Also suitable for large animals such as pigs, sheeps or horses, dogs, cats or rabbits

Type	Dimensions (L x W x H) in mm	Order no.
FeedMe roll	~ 270 x 110 x 110	3084300
FeedMe Futterrolle reinforced	~ 270 x 110 x 110	3084301

BRACKETS FOR FEEDME

Various mounting brackets have been developed for permanent installation on animal enclosures. Depending on the type, the feed rolls can be attached to the outside or inside of the enclosure using vertical or horizontal bars, or even to a wire mesh panel. This makes the feed roll a useful enrichment tool even for large animals such as horses, sheep, pigs or primates.



Type	Dimensions (L x W x H) in mm	Order no.
FeedMe roll bracket - NHP Kit (reinforced)		3084302
FeedMe roll bracket - Pig Kit (reinforced)		3084303
FeedMe bracket for attaching to the fence		3084304