

GLOVES BITE RESISTANT



Order Number on request

Protect your fingers from mice, rats hamsters or Guinea pigs teeth

In order to protect the thumb, index finger and middle finger that are easily bitten by mice and rats, are coated with the material KEVLAR.

- Bite and cut prevention**

Size	S	M
	20 pairs / bag	20 pairs / bag

When young scientists start working with small rodents fears of beeing bitten can arise. Here, the bite-resistant gloves help to provide the same tactile sensation during handling as when using bare hands. In this way, the bite-resistant gloves can help to minimise stress for both humans and animals. For completely hygienic working conditions, simply pull standard latex gloves over the top.

- The cut resistance of this product means that it is difficult to cut. It does not mean puncture resistance.
- KEVLAR is discolored by sunlight and UV light. If the color changes, the strength of the material may be reduced. It is therefore recommended to store it within the delivered aluminium bag in which the product is delivered.
- Do not bleach when washing. After washing dry in the shadow without exposing to direct sunlight.